

INSTRUCTIONS FOR USE

Disposable Medical Face Mask

【Product Name】

Disposable Medical Face Mask

【Brief Introduction】

Melt blown non-woven fabric can filter harmful substances in the air and effectively control the air of harmful substances to be filtered through the filter layer of the mask before being inhaled or exhaled by person.

【Model description】

70010001, 70011001

【Intended Use】


It is intended to cover the user's mouth, nose and jaw. It is used to wear and block the disposable mask that exhales or ejects pollutants from the mouth and nose in the general medical environment.

【Contraindications】

Use it with caution if you are allergic to non-woven fabrics.

【Structure composition】

Disposable medical face mask is composed of a face mask body, a nose clip, and the face mask body has three layers of structure. The structure diagram is as follows:



1: Mask body 2: A nose clip 3: A mask belt

Fig 1 Flat earhook type

【Raw materials】

The inner and outer layers are polypropylene spunbond non-woven fabric, and the middle layer is composed of melt blown non-woven fabric, plastic nose clip and elastic mask belt

【Sterilization state】



The product is non-sterile.

【Precautions and warnings】

- 1) Please read the instructions carefully to avoid touching the inner layer of the mask;
- 2) This product has not been sterilized, please check the packaging before use, if the packaging is damaged, it is forbidden to use;
- 3) This product is disposable and cannot be reused;
- 4) In case of liquid splashing, mask damage, damp and obvious increase of respiratory resistance during wearing, please replace it in time;
- 5) After the product is used, it shall be handled according to the requirements of the hospital or environmental protection department.

【Steps for usage】

- 1) Open the package, take out the mask, hold both ends of the mask, with the nose bar facing up;
- 2) The mask covers the mouth, nose and mandible, and the ear belt is placed behind the ear;
- 3) Adjust the bridge of the nose to fit it with the bridge of the nose;
- 4) Properly adjust the position of the mask to make it comfortable to wear;
- 5) Wear the mask correctly according to the front and back signs. Do not touch the inside of the mask with your hands.

【Performance requirements】

The performance complies with EN 14683:2019+AC: 2019. BFE \geq 95%

【Product List】

No	Name
1	Disposable Medical Face Mask
2	Manual

【Product storage and transportation】

- 1) During transportation, the product shall be protected from heavy pressure, direct sunlight and rain and snow.
- 2) Store in a ventilated, dry, non-corrosive gas environment. Keep away from fire and flammable materials.

【Expiration date】

The product is valid for 2 years from the production date, please use it within the validity period.

【Definitions of Signs】



	" CE " MARKING		CAUTION
	DO NOT USE IF PACKAGE IS DAMAGED		DO NOT REUSE
	KEEP AWAY FROM SUNLIGHT		KEEP DRY
	BATCH CODE		USE BY
	DATE OF MANUFACTURE		CONSULT INSTRUCTIONS FOR USE
	MANUFACTURER		EUROPEAN REPRESENTATIVE
	Non-sterile		Medical device

**Manufacturer**

Name: Shanghai Oma Safety Products Co., Ltd.

ADD: Room 506, No. 23, Lane 423, Xincun Road, Shanghai 200065, China

Zip Code: 200065

Tel: +86 (21)-32092351

Fax: +86-21-32092350

E-mail: sales-1@oma-safety.cn

**Authorized representative in the European Community**

Name: CMC MEDICAL DEVICES & DRUGS S.L.

Add: C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain

E-mail: victoriali@cmcmmedicaldevices.com

【Version】 2020.05.10**【Revision No】** A/0